



Instructions

FOR IN VITRO DIAGNOSTIC
USE ONLY

INTENDED USE

The DrugSmart Test Cup is a one-step immunoassay for the qualitative detection of multiple drugs and drug metabolites in human urine at the following cutoff concentrations:

| Test | Calibrator | Cut-off (ng/ml) |
|---------|--------------------------------|-----------------|
| AMP | Amphetamine | 1000 |
| BAR | Secobarbital | 300 |
| BUP | Buprenorphine | 10 |
| BZO | Oxazepam | 300 |
| COC150 | Benzoylcegonine | 150 |
| COC300 | Benzoylcegonine | 300 |
| MDMA | 3,4-methylenedioxyamphetamine | 500 |
| MET500 | Methamphetamine | 500 |
| MET1000 | Methamphetamine | 1000 |
| MTD | Methadone | 300 |
| OPI300 | Morphine | 300 |
| OPI2000 | Morphine | 2000 |
| OXY | Oxycodone | 100 |
| PCP | Phencyclidine | 25 |
| PPX | Propoxyphene | 300 |
| TCA | Nortriptyline | 1000 |
| THC | 11-nor- Δ^9 -THC-9-COOH | 50 |

The configurations of the DrugSmart Test Cup consist of any combination of the drugs listed above with or without specimen validity test. The specimen validity test provides information regarding the integrity of urine sample in the drugs of abuse test by the semi-quantitative determination of pH, Specific Gravity and Oxidants in human urine. The DrugSmart Test Cup is used to obtain a visual, qualitative result and is intended for professional use only.

This assay provides only a preliminary result. Clinical consideration and professional judgment must be applied to any drug of abuse test result, particularly in evaluating a preliminary positive result. In order to obtain a confirmed analytical result, a more specific alternate chemical method is needed. Gas Chromatography/Mass Spectroscopy (GC/MS) is the preferred confirmation method.

SUMMARY AND EXPLANATION

Amphetamine/Methamphetamine, amphetamine, and metabolites are potent central nervous system stimulants. Acute higher doses induce euphoria, alertness, and sense of increased energy and power. More acute responses produce anxiety, paranoia, psychotic behavior, and cardiac dysrhythmias. Methamphetamine is excreted in urine as amphetamine and oxidized as deaminated and hydroxylated derivatives. However, methamphetamine is also excreted to some extent unchanged. Thus the presence of the parent compound in the urine indicates methamphetamine use.

Barbiturates are classified as central nervous system depressants. These products produce a state of intoxication that is similar to alcohol intoxication. Symptoms include slurred speech, loss of motor coordination and impaired judgment. Depending on the dose, frequency, and duration of use, one can rapidly develop tolerance, physical dependence and psychological dependence on barbiturates. Barbiturates are taken orally, or by intravenous and intramuscular injections. They are excreted in urine as parent compound as well as metabolites.

Benzodiazepines are central nervous system (CNS) depressants commonly prescribed for the short-term treatment of anxiety and insomnia. In general, benzodiazepines act as hypnotics in high doses, as anxiolytics in moderate doses and as sedatives in low doses. The use of benzodiazepines can result in drowsiness and confusion. Psychological and physical dependence on benzodiazepines can develop if high doses of the drug are given over a prolonged period. Benzodiazepines are taken orally or by intramuscular or intravenous injection, and are extensively oxidized in the liver to metabolites. Parent compounds, as well as metabolites are excreted in the urine.

Buprenorphine is a synthetic thebaine derivative that has both analgesic and opioid antagonist properties. As an analgesic, it is about 25 to 40 times more potent than morphine. Symptoms of overdosage include confusion, dizziness, pinpoint pupils, hallucinations, hypotension, respiratory difficulty, seizures and coma. Buprenorphine is metabolized in man primarily by N-dealkylation and conjugation to form norbuprenorphine, which is pharmacologically active, and conjugates of buprenorphine and norbuprenorphine. Within 144 hours of a single intramuscular dose of drug, 95% is eliminated, with 68% in the feces and 27% in the urine. Buprenorphine and its metabolites in urine may be detected as a result of buprenorphine, norbuprenorphine, Buprenorphine-3-beta-D-glucuronide, and Norbuprenorphine-3-beta-D-glucuronide.

Cocaine is a potent central nervous system stimulant and a local anesthetic found in the leaves of the coca plant. The psychological effects induced by using

cocaine are euphoria, confidence and sense of increased energy. These psychological effects are accompanied by increased heart rate, dilation of the pupils, fever, tremors and sweating. Cocaine is excreted in the urine primarily as benzoylecgonine in a short period of time. Benzoylecgonine has a biological half-life of 5 to 8 hours, which is much longer than that of cocaine (0.5 to 1.5 hour), and can be generally detected for 24 to 60 hours after cocaine use or exposure.

3,4-methylenedioxyamphetamine is classified as both a stimulant and a hallucinogen. Like methamphetamine, adverse effects of 3,4-methylenedioxyamphetamine use include jaw clenching, teeth grinding, dilated pupils, perspiring, anxiety, blurred vision, vomiting, and increased blood pressure and heart rate. Overdose of 3,4-methylenedioxyamphetamine may cause heart failure or extreme heart stroke. 3,4-methylenedioxyamphetamine is taken orally in tablets or capsules and excreted in urine as parent compound as well as metabolic.

Methadone is a synthetic analgesic drug originally used for the treatment of narcotic addiction. The psychological effects induced by using methadone are analgesia, sedation, and respiratory depression. Overdose of methadone may cause coma or even death. Methadone is taken orally or intravenously and is metabolized in the liver and has a biological half-life of 15-60 hours.

Opiates, such as heroin, morphine, and codeine, are central nervous system (CNS) depressants. The use of opiates at high doses produces euphoria and release from anxiety. Physical dependence is apparent in users and leads to depressed coordination, disrupted decision making, decreased respiration, hypothermia and coma. Heroin is quickly metabolized to morphine, morphine glucuronide and 6-acetylmorphine. Thus, the presence of morphine (or the metabolite, morphine glucuronide) in the urine indicates heroin, morphine, and/or codeine use.

Oxycodone is a semi-synthetic opioid with a structural similarity to codeine. It produces potent euphoria, analgesic and sedative effects, and has a dependence liability similar to morphine. Oxycodone is most often administered orally and is metabolized by demethylation to noroxycodone and oxymorphone followed by glucuronidation and excreted in urine. The window of detection for oxycodone in urine is expected to be similar to that of other opioids such as morphine.

Phencyclidine, commonly known as "angel dust" and "crystal cyclone", is an arylcyclohexylamine that is originally used as an anesthetic agent and a veterinary tranquilizer. The drug is abused by oral or nasal ingestion, smoking, or intravenous injection. It produces hallucinations, lethargy, disorientation, loss of coordination, trance-like ecstatic states, a sense of euphoria and visual distortions. It is well absorbed following all routes of administration. Unchanged PCP is excreted in urine in moderate amounts (10% of the dose).

Propoxyphene is a mildly effective narcotic analgesic that has been in clinical use. It is less potent than codeine, and bears a close structural relationship to methadone. Propoxyphene is available in oral formulations either as the hydrochloride (32 or 65 mg) or as the napsylate salt (50 or 100 mg), and is often found in combination with aspirin or acetaminophen. Overdosage with propoxyphene can result in stupor, coma, convulsions, respiratory depression, cardiac arrhythmias, hypotension, pulmonary edema and circulatory collapse. Propoxyphene is metabolized primarily via N-demethylation to norpropoxyphene. The amounts of metabolites excreted in the 20 hour urine following a 130 mg single oral dose of propoxyphene hydrochloride were: 1.1% propoxyphene, 13.2% norpropoxyphene and 0.7% dinorpropoxyphene.

Tetrahydrocannabinol is generally accepted to be the principle active component in marijuana. When ingested or smoked, it produces euphoric effects. Abusers exhibit central nervous system effects, altered mood and sensory perceptions, loss of coordination, impaired short term memory, anxiety, paranoia, depression, confusion, hallucinations and increased heart rate. When marijuana is ingested, the drug is metabolized by the liver, the primary metabolite of marijuana excreted in the urine is 11-nor- Δ^9 -tetrahydrocannabinol-9-carboxylic acid. Therefore, the presence of detected cannabinoids, including the primary carboxyl metabolite, in the urine indicate marijuana/cannabis use.

Tricyclic antidepressants (TCAs) have been prescribed for depression and compulsive disorders. Because of the possibility of causing serious cardiac complications, TCAs can be lethal if misused at high doses. TCAs are taken orally or sometimes by injection. TCAs are metabolized in the liver. Both TCAs and their metabolites are excreted in urine mostly in the form of metabolites for up to ten days.

The length of time following drug use of which a positive result may occur is dependent upon several factors, including the frequency and amount of drug, metabolic rate, excretion rate, drug half-life, and the drug user's age, weight, activity and diet.

Adulteration of urine samples may cause erroneous results in drugs of abuse test by either interfering with the drug screening test and/or destroying the drugs in the urine. Dilution of urine with water is probably the simplest urine adulteration method. Bleach, vinegar, urine, sodium bicarbonate, sodium nitrite, Drano, soft drinks and hydrogen peroxide are the examples of adulterants used to adulterate the urine sample. It is important to insure the integrity of urine samples in drugs of abuse test.

TEST PRINCIPLE

The DrugSmart Test Cup is based on the principle of competitive immunochemical reaction between a chemically labeled drug (drug-protein conjugate) and the drug or drug metabolites, which may be present in the urine sample for the limited antibody binding sites. The test contains

a nitrocellulose membrane strip pre-coated with drug-protein conjugate in the test region and a pad containing colored antibody-colloidal gold conjugate. During the test, the urine sample is allowed to migrate upward and dehydrate the antibody-colloidal gold conjugate. The mixture then migrates along the membrane chromatographically by the capillary action to the immobilized drug-protein band on the test region. When drug is absent in the urine, the colored antibody-colloidal gold conjugate and immobilized drug-protein bind specifically to form a visible line in the test region as the antibody complexes with the drug-protein. When drug is present in the urine, it will compete with drug-protein for the limited antibody sites. The line on the test region will become less intense with increasing drug concentration. When a sufficient concentration of drug is present in the urine, it will fill the limited antibody binding sites. This will prevent attachment of the colored antibody-colloidal gold conjugate to the drug-protein on the test region. Therefore, the presence of the line on the test region indicates a **negative** result for the drug and the absence of the test line on the test region indicates a **positive** result for the drug.

A visible line generated by a different antigen/antibody reaction is also present at the control region of the test strip. This line should always appear, regardless of the presence of drugs or metabolites in the urine sample. This means that a **negative** urine sample will produce both test line and control line, and a **positive** urine sample will generate only control line. The presence of control line serves as a built-in control, which demonstrates that the test is performed properly.

The DrugSmart Test Cup for the specimen validity test is based on the color response of chemical indicators at presence of adulterants. Specific Gravity, pH and Oxidants (OX) are tested to determine the integrity of urine samples.

pH: The pH determination of urine sample is based on color change of indicator in different acidic or basic medium. The normal urine pH ranges from 4.5 to 8. Urine pH below 4 or above 9 indicates adulteration acidic.

OX: Bleach or other oxidizing agents react with an oxidant indicator to form a color complex. Observation of a blue-green, brown, or orange color indicates adulteration with bleach or other oxidizing agents.

S.G.: Random urine may vary in specific gravity from 1.003 – 1.030. Normal adults with normal diets and normal fluid intake will have an average urine specific gravity of 1.016 – 1.022. Elevated urine specific gravity values may be obtained in the presence of moderate quantities of protein. A urine specimen with a specific gravity level of less than 1.003 can be an indication of substitution.

REAGENTS & MATERIALS SUPPLIED

- 25 individually wrapped test cups. Each device consists of different test strips in a plastic test strip holder. The test strip contains a colloidal gold pad coated with antibody (or drug-protein) and rabbit antibody. It also contains a membrane coated with drug-bovine protein conjugate (or antibody) in the test band and goat anti-rabbit antibody in the control band and adulterant pads when applicable.
- One instruction sheet.
- 25 Security Seals.
- One Adulteration Color Comparison Chart for interpretation of adulteration test result (when applicable).

MATERIAL REQUIRED BUT NOT PROVIDED

- Timer

WARNINGS AND PRECAUTIONS

- For professional *in vitro* diagnostic use only
- Urine specimens may be potentially infectious. Proper handling and disposal methods should be established.
- Avoid cross-contamination of urine samples by using a new specimen collection container for each urine sample.
- Test device should remain sealed until ready for use.
- Do not use the test kit after the expiration date.
- A positive test result does not always mean an individual has taken the drug illegally as the drug may have been administered legally.

STORAGE

The DrugSmart Test Cup should be stored at 2-30°C (36-86°F) in the original sealed pouch. Do not freeze. Do not store and or expose reagent kits at temperature greater than 30°C.

SPECIMEN COLLECTION AND HANDLING

Fresh urine does not require any special handling or pretreatment. A fresh urine sample should be collected in the test cup with a minimum of 30ml volume. **The DrugSmart Test Cup employs a thermal strip to validate the urine collection. This device should be checked immediately after collection.**

Note: Urine specimens can be transferred from a urine collection container into the DrugSmart Test Cup, if necessary.

TEST PROCEDURE

If specimen, control, or test devices have been stored at refrigerated temperatures, allow them to warm to room temperature before testing.

Do not open test device pouch until ready to perform the test.


1. Tear open the foil pouch and remove the Test Cup.
2. Issue the device to the individual to be tested.
3. Have them urinate directly into the Test Cup. Ensure the specimen is above the minimum level line indicated on the test cup label.
4. The cup must be returned immediately to the collector. Authorized personnel at the collection site is to remove tear-off label and read the results at five minutes post collection.

NOTE: In order to prevent any incorrect results, the test results should **not** be interpreted after 10 minutes.

If the test cup includes an adulteration panel, read the results of adulteration test immediately, prior to the drug test result, by visually comparing the color of reagent pads to the corresponding blocks on the Color Chart at the time indicated.

If results show that the sample was adulterated, DO NOT read the drug test results.

INTERPRETATION OF RESULTS



Negative (-): Colored lines appear in both Control Region (C) and Test Region (T). The line in the control region is the control line, which is used to indicate proper performance of the device. The line in the test region is the drug probe line. The test line may have varying intensity either weaker or stronger in color than that of the control line. A negative result for a drug indicates that the concentration of that drug in urine is below the cutoff level.

Positive (+): Colored line appears in the control region only. No line appears in the test region. The complete absence of a test line indicates a positive result for that drug. A preliminary positive result for a drug indicates that the concentration of that drug in urine is at or above the cutoff level.

Invalid: No colored line appears in the control region. If the control line does not form, the test result is inconclusive and should be repeated. Each strip in the DrugSmart Test Cup is read and functions independently. An invalid result on one test strip does not invalidate other results derived from the same device.

Each test strip is read individually and independently of one another.

Specimen Validity Tests:

Specimen validity test results are obtained by directly comparing the color of each test pad with the color block of Adulteration Color Comparison Chart. Adulterated urine sample will produce abnormal color response. Unadulterated urine sample will produce normal color response.

QUALITY CONTROL

An internal procedural control is included in the test device. A line must form in the Control band region regardless of the presence or absence of drugs or metabolites. The presence of the line in the Control region indicates that the proper sample volume has been used and that the reagents are migrating properly. If the line in the Control region does not form, the test is considered invalid.

To ensure proper kit performance, it is recommended that the test devices be tested once a week with external controls. External controls are available from commercial sources. It is important to make sure that the control values are within established limits. If the values of external control do not fall within established limits, the test results are invalid. Additional controls may be tested according to guidelines or requirements of local, state, and/or federal regulations or accrediting organizations.

LIMITATIONS OF PROCEDURE

- The assay is designed for use with human urine only.
- A positive result with any of the tests indicates only the presence of a drug/metabolite and does not indicate or measure intoxication.
- There is a possibility that technical or procedural error as well other substances as factors not listed may interfere with the test and cause false results. See SPECIFICITY for lists of substances that will produce positive results, or that do not interfere with test performance.
- If adulteration is suspected, the test should be repeated with new sample.

PERFORMANCE CHARACTERISTICS

Accuracy

The accuracy of the DrugSmart Test Cup was evaluated in comparison to commercially available drug screen tests. Sixty (60) negative urine samples collected from presumed non-user volunteers were tested by both DrugSmart Test Cup and commercially available drug screen tests. Of these negative urine samples tested, all were found negatives by both methods. In a separate study, positive urine samples, obtained from clinical laboratories where the drug concentrations were determined by GC/MS (TCA concentrations were determined by HPLC), were tested by DrugSmart Test Cup and commercial drug screen tests. The results of accuracy study are presented in the following table:

| Drug Test | | GC/MS (<-50% C/O) | GC/MS (-50% C/O to C/O) | GC/MS (C/O to +50% C/O) | GC/MS (>+50% C/O) | % Agreement with GC/MS |
|-----------|-----|----------------------|-------------------------------|-------------------------------|----------------------|---------------------------|
| AMP | (+) | 0 | 0 | 10 | 55 | 98.5 |
| | (-) | 15 | 9 | 1 | 0 | 100 |
| BAR | (+) | 0 | 1 | 5 | 83 | 97.8 |
| | (-) | 15 | 7 | 2 | 0 | 95.7 |
| BUP | (+) | 0 | 0 | 8 | 35 | 97.7 |
| | (-) | 18 | 6 | 1 | 0 | 100 |
| BZO | (+) | 0 | 2 | 13 | 37 | 100 |
| | (-) | 18 | 18 | 0 | 0 | 94.7 |
| COC150 | (+) | 0 | 1 | 7 | 60 | 100 |
| | (-) | 15 | 10 | 0 | 0 | 96.2 |
| COC300 | (+) | 0 | 0 | 8 | 71 | 98.8 |
| | (-) | 15 | 8 | 1 | 0 | 100 |
| MDMA | (+) | 0 | 1 | 6 | 37 | 100 |
| | (-) | 24 | 6 | 0 | 0 | 96.8 |
| MET500 | (+) | 0 | 2 | 8 | 64 | 100 |
| | (-) | 15 | 4 | 0 | 0 | 90.5 |
| MET1000 | (+) | 0 | 0 | 5 | 58 | 98.4 |
| | (-) | 20 | 8 | 1 | 0 | 100 |
| MTD | (+) | 0 | 0 | 6 | 65 | 98.6 |
| | (-) | 15 | 5 | 1 | 0 | 100 |
| OPI300 | (+) | 0 | 1 | 6 | 77 | 100 |
| | (-) | 16 | 6 | 0 | 0 | 95.7 |
| OPI2000 | (+) | 0 | 2 | 9 | 45 | 100 |
| | (-) | 15 | 6 | 0 | 0 | 91.3 |
| OXY | (+) | 0 | 1 | 6 | 47 | 100 |
| | (-) | 15 | 7 | 0 | 0 | 95.7 |
| PCP | (+) | 0 | 0 | 4 | 56 | 96.8 |
| | (-) | 15 | 4 | 2 | 0 | 100 |

| | | | | | | |
|-----|-----|----|----|----|----|------|
| PPX | (+) | 0 | 0 | 6 | 64 | 98.6 |
| | (-) | 10 | 7 | 1 | 0 | 100 |
| TCA | (+) | 0 | 1 | 12 | 9 | 100 |
| | (-) | 23 | 11 | 0 | 0 | 97.1 |
| THC | (+) | 0 | 1 | 24 | 32 | 100 |
| | (-) | 15 | 12 | 0 | 0 | 96.4 |

Precision

The precision of the DrugSmart Test Cup was evaluated by testing three lots of the test devices at four study sites with spiked drug sample solutions on three consecutive days. Sample concentrations were confirmed by GC/MS.

| | | | | | | |
|-----------------|-------|-------|--------|-------|--------|-------|
| AMP (ng/ml) | 0 | 500 | 750 | 1000 | 1250 | 1500 |
| (+/-) | 0/135 | 0/135 | 34/101 | 75/60 | 110/25 | 135/0 |
| BAR (ng/ml) | 0 | 150 | 225 | 300 | 375 | 450 |
| (+/-) | 0/135 | 0/135 | 34/101 | 74/61 | 102/33 | 135/0 |
| BUP (ng/ml) | 0 | 5 | 7.5 | 10 | 12.5 | 15 |
| (+/-) | 0/135 | 0/135 | 33/102 | 73/61 | 101/34 | 135/0 |
| BZO (ng/ml) | 0 | 150 | 225 | 300 | 375 | 450 |
| (+/-) | 0/135 | 0/135 | 29/106 | 75/60 | 107/28 | 135/0 |
| COC 150 (ng/ml) | 0 | 75 | 112.5 | 150 | 187.5 | 225 |
| (+/-) | 0/135 | 0/135 | 33/102 | 75/60 | 105/30 | 135/0 |
| COC300 (ng/ml) | 0 | 150 | 225 | 300 | 375 | 450 |
| (+/-) | 0/135 | 0/135 | 30/105 | 65/70 | 96/36 | 135/0 |
| MDMA (ng/ml) | 0 | 250 | 375 | 500 | 625 | 750 |
| (+/-) | 0/135 | 0/135 | 35/100 | 75/60 | 95/40 | 135/0 |
| MET500 (ng/ml) | 0 | 250 | 375 | 500 | 625 | 750 |
| (+/-) | 0/135 | 0/135 | 32/103 | 77/58 | 99/36 | 135/0 |
| MET1000 (ng/ml) | 0 | 500 | 750 | 1000 | 1250 | 1500 |
| (+/-) | 0/135 | 0/135 | 31/104 | 77/58 | 98/37 | 135/0 |
| MTD (ng/ml) | 0 | 150 | 225 | 300 | 375 | 450 |
| (+/-) | 0/135 | 0/135 | 31/104 | 69/66 | 95/40 | 135/0 |
| OPI300 (ng/ml) | 0 | 150 | 225 | 300 | 375 | 450 |
| (+/-) | 0/135 | 0/135 | 33/102 | 70/65 | 95/40 | 135/0 |
| OPI2000 (ng/ml) | 0 | 1000 | 1500 | 2000 | 2500 | 3000 |
| (+/-) | 0/135 | 0/135 | 37/98 | 76/59 | 104/31 | 135/0 |
| OXY (ng/ml) | 0 | 50 | 75 | 100 | 125 | 150 |
| (+/-) | 0/135 | 0/135 | 50/85 | 86/49 | 111/24 | 135/0 |
| PCP (ng/ml) | 0 | 12.5 | 18.75 | 25 | 31.25 | 37.5 |
| (+/-) | 0/135 | 0/135 | 26/109 | 62/73 | 99/36 | 135/0 |
| PPX (ng/ml) | 0 | 150 | 225 | 300 | 375 | 450 |
| (+/-) | 0/135 | 0/135 | 34/101 | 77/58 | 103/32 | 135/0 |
| TCA (ng/ml) | 0 | 500 | 750 | 1000 | 1250 | 1500 |
| (+/-) | 0/135 | 0/135 | 24/111 | 60/75 | 99/36 | 135/0 |
| THC (ng/ml) | 0 | 25 | 37.5 | 50 | 62.5 | 75 |
| (+/-) | 0/135 | 0/135 | 27/108 | 58/77 | 91/44 | 135/0 |

Specificity

The specificity for the DrugSmart Test Cup was determined by testing various drugs, drug metabolites, and other compounds that are likely to be present in urine. All compounds were prepared in drug-free normal human urine.

The following compounds produce positive results when tested at levels greater than the concentrations listed below.

| Compound | Conc. (ng/ml) | Compound | Conc. (ng/ml) |
|---------------------|---------------|-------------------|---------------|
| Amphetamine | | | |
| d-Amphetamine | 1,000 | d-Methamphetamine | 50,000 |
| dl-Amphetamine | 2,500 | (+/-)3,4-MDMA | 50,000 |
| (+/-)3,4-MDA | 1,250 | | |
| Barbiturates | | | |
| Secobarbital | 300 | Butobarbital | 400 |
| Allobarbital | 600 | Butalbital | 300 |
| Alphenal | 200 | Butethal | 450 |
| Amobarbital | 1500 | Pentobarbital | 400 |
| Aprobarbital | 300 | Phenobarbital | 450 |
| Barbital | 1500 | | |

| | | | |
|----------------------------------|----------|---|----------|
| Benzodiazepines | | | |
| Oxazepam | 300 | Flunitrazepam | 300 |
| Alprazolam | 400 | Flurazepam | 300 |
| Bromazepam | 250 | Lorazepam | 500 |
| Chlordiazepoxide | 300 | Medazepam | 300 |
| Clobazam | 1000 | Nitrazepam | 250 |
| Clonazepam | 500 | Nordiazepam | 150 |
| Clorazepate | 150 | Prazepam | 500 |
| Desalkylflurazepam | 200 | Temazepam | 200 |
| Diazepam | 450 | Triazolam | 450 |
| Estazolam | 300 | | |
| Buprenorphine | | | |
| Buprenorphine | 10 | Buprenorphine-3- beta-D-glucuronide | 7.5 |
| Norbuprenorphine | 25,000 | | |
| Codeine | >100,000 | | |
| Morphine | >100,000 | Norbuprenorphine-3- beta-D-glucuronide | 150 |
| Nalorphine | 10,000 | | |
| Cocaine (150) | | | |
| Benzoylcegonine | 150 | Cocaethylene | >100,000 |
| Cocaine | 5,000 | Ecgonine methyl esters | >100,000 |
| Ecgonine | >100,000 | | |
| Cocaine (300) | | | |
| Benzoylcegonine | 300 | Cocaine | 300 |
| Methamphetamine (500) | | | |
| d-Methamphetamine | 500 | (+/-)3,4-MDMA | 2,000 |
| d-Amphetamine | 50,000 | l-Methamphetamine | 10,000 |
| l-Amphetamine | >100,000 | Ephedrine | 100,000 |
| (+/-)3,4-MDEA | 50,000 | Mephentermine | 50,000 |
| (+/-)3,4-MDA | 100,000 | | |
| Methamphetamine (1000) | | | |
| d-Methamphetamine | 1000 | (+/-)3,4-MDMA | 3,000 |
| d-Amphetamine | 50,000 | l-Methamphetamine | 10,000 |
| l-Amphetamine | >100,000 | Ephedrine | >100,000 |
| (+/-)3,4-MDEA | 50,000 | Mephentermine | 75,000 |
| (+/-)3,4-MDA | 100,000 | | |
| MDMA | | | |
| (+/-)3,4-MDMA | 500 | (+/-)3,4-MDA | 4,000 |
| (+/-)3,4-MDEA | 450 | | |
| Methadone | | | |
| (+/-) Methadone | 300 | Methadol | 1,500 |
| Opiates (300) | | | |
| Morphine | 300 | Hydrocodone | 500 |
| Codeine | 300 | Hydromorphone | 500 |
| Ethylmorphine | 300 | Morphine-3-glucuronide | 300 |
| Heroin | 750 | Nalorphine | 5,000 |
| Opiates (2000) | | | |
| Morphine | 2,000 | Hydrocodone | 4,000 |
| Codeine | 2,000 | Hydromorphone | 5,000 |
| Ethylmorphine | 1,000 | Morphine-3-glucuronide | 2,500 |
| Heroin | 5,000 | Nalorphine | 5,000 |
| Oxycodone | | | |
| Oxycodone | 100 | Morphine | >100,000 |
| Hydrocodone | 5000 | Codeine | 50,000 |
| Hydromorphone | 50,000 | Heroin | >100,000 |
| PCP | | | |
| Phencyclidine | 25 | Tenocyclidine | 2,000 |
| PPX | | | |
| D-Propoxyphene | 300 | D-Norpropoxyphene | 300 |
| THC | | | |
| 11-nor- Δ^9 -THC-9-COOH | 50 | Δ^9 -tetrahydrocannabinol | 5,000 |
| 11-hydroxy- Δ^9 -THC | 1,000 | Cannabinol | 10,000 |
| Δ^8 -tetrahydrocannabinol | 5,000 | Cannabidiol | >100,000 |

Tricyclic Antidepressant

| | | | |
|---------------|-------|-------------|-------|
| Nortriptyline | 1,000 | Promazine | 1,500 |
| Nordoxepin | 2,000 | Desipramine | 400 |
| Trimipramine | 2,000 | Doxepin | 3,000 |
| Amitriptyline | 1,500 | Maprotiline | 2,000 |

Interference

Two pools of drug-free urine were spiked with drug standards to 50% below and 50% above cutoff concentrations. The drug concentrations were confirmed by GC/MS. The following compounds were evaluated for potential positive and/or negative interference with the DrugSmart Test Cup. All compounds were dissolved in the spiked sample solutions and tested with DrugSmart Test Cup. An unaltered sample was used as a control. No positive interference or negative interference was found for the following compounds when tested at concentrations up to 100 µg/ml.

| | | |
|------------------------|-------------------------|--------------------------------|
| Acetaminophen | Diphenhydramine | (+/-)-Norephedrine |
| Acetone | Dopamine | Oxalic Acid |
| Albumin | (+/-)-Epinephrine | Penicillin-G |
| Ampicillin | Erythromycin | Pheniramine |
| Ascorbic Acid | Ethanol | Phenothiazine |
| Aspartame | Furosemide | l-Phenylephrine |
| Aspirin | Glucose | β -Phenylethylamine |
| Atropine | Guaiacol Glyceryl Ether | Procaine |
| Benzocaine | Hemoglobin | Quinidine |
| Bilirubin | Ibuprofen | Ranitidine |
| Caffeine | (+/-)-Isoproterenol | Riboflavin |
| Chloroquine | Ketamine | Sodium Chloride |
| (+)-Chlorpheniramine | Levorphanol | Sulindac |
| (+/-)-Chlorpheniramine | Lidocaine | Theophylline |
| Creatine | (+)-Naproxen | Tyramine |
| Dexbrompheniramine | Niacinamide | 4-Dimethylaminoantipyrine |
| Dextromethorphan | Nicotine | (1R,2S)-(-)-N-Methyl-Ephedrine |

Effect of Specimen pH

Drug sample solutions with 50% below and 50% above cutoff concentrations were adjusted to pH 4-9 and tested using DrugSmart Test Cup. An unaltered sample was used as a control. The results demonstrate that varying ranges of specimen pH do not interfere with the performance of the test.

Effect of Specimen Specific Gravity

Drug sample solutions with 50% below and 50% above cutoff concentrations were adjusted to specific gravity 1.003-1.04 and tested using DrugSmart Test Cup. An unaltered sample was used as a control. The results demonstrate that varying ranges of specimen specific gravity do not interfere with the performance of the test.

BIBLIOGRAPHY OF SUGGESTED READING

- Baselt, R. C., Disposition of Toxic Drugs and Chemicals in Man, Biomedical Publications, Davis, CA, 1982.
- Urine testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA), Research Monograph 73, 1986.
- Fed. Register, Department of Health and Human Services, Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970-11979, 1988.
- Liu, Ray H. and Goldberger, Bruce A., Handbook of Workplace Drug Testing, AACCC Press (1995).
- Gilman, A. G. and Goodman, L. S., The Pharmacological Basis of Therapeutics, eds. MacMillan Publishing, New York, NY, 1980.